Compact

VBH CE

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
COMPANY
WITH QUALITY SYSTEM
CERTIFIED BY DNV
= ISO 9001/2000 =

MICROBIOLOGICAL SAFETY CABINET

INTRODUCTION

VBH Compact are vertical laminar flow biohazard cabinets, Class II A.

VBH Compact are designed for the situation where is requested the operator, product and environmental protection from dangerous effects due to uncontrolled diffusion of air-transported contaminants and, in the same time, to avoid any contamination from the environment to the product during its handling.

STANDARDS AND CERTIFICATIONS

VBH Compact cabinets produced by Steril are designed to fulfil and satisfy the European Norm EN12469 and country standards (presently superseded by EN12469) like German DIN 12950/10, French NF44-201 and English BS 5726.

FUNCTIONING

VBH Compact cabinets are manufactured with a frontal working opening for the aspiration of 30% of the total air involved, at a minimum velocity of 0,40 m/s, and the remaining 70% is recirculated through a HEPA filter in the working chamber with a mean laminar flow velocity of 0,40 m/s (as defined by European standard EN 12469: 2000).

The containment index, evaluated on the total surface of the front working opening, corresponds to APF (Aperture Protection Factor) equal or more than 1×10^5 according EN12469.

The microprocessor monitors in real time all critical parameters of the unit (laminar air flow velocity, exhaust flow rate and position of the front screen) and controls the optimal function of the fans even if in presence of effects caused by progressive clogging of the HEPA filters.

TECHNICAL SPECIFICATION

Work chamber and spillage tray in stainless steel AISI 304L * 2B* finishing.

Working surface in stainless steel AISI 304L *2B* finishing, removable and autoclavable in perforated and solid version.

Front window, in laminated safety glass electrical driven. Window is desing in order to obtain an aerosol-light in two positions (working position and close position, with automatic switch off of the motor blower).

Ergonomic front closure panel in anodised aluminium complete of germicidal UV lamp (optional).

Lighting on working surface by means of 2 fluorescent electronic lamps fitted in non-contaminated area.

Construction structure in steel epoxy-polyester powder coated. **Ventilation system** by means of n°1 IP55 protection factor motor blower.

Recirculating and exhaust **air filtration** achieved by **H14 HEPA filters as per EN 1822 (with efficiency 99.995% MPPS)**. Maintenance completely from the front of the cabinet.

SAFETY DEVICES

Microprocessor monitoring with alarms of the following parameters:

- Exhaust flow rate out of the pre-set range
- Laminar air flow velocity out of the pre-set range
- Front window in wrong position
- Fans failure.

STANDARD FEATURES

- Decontamination program
- UV timer program
- Power socket 4A IP55
- Vacuum air line
- Fuel gas line
- Multilanguage display

OPTIONAL

- UV germicidal lamp
- PC interface for monitoring and remote servicing with lithium battery
- Double HEPA exhaust filter
- Activated carbon exhaust filter
- Support stand
- Chest of drawers on wheels

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